

LX



Gips-Instrumente  
Plaster Instruments

## Gips-Instrumente – Plaster Instruments

Gips-Instrumente  
Plaster Instruments

LX



1/2

ESMARCH

LX505R

Gipsmesser  
Plaster knives  
180 mm, 7"



1/2

ENGEL

LX522C

Gipssäge  
Plaster saw  
150 mm, 6"



1/2

LISTER

BC885R-BC977R

BC885R	schwarz	black	BC973R	gelb	yellow
BC970R	rot	red	BC975R	orange	orange
BC971R	blau	blue	BC977R	pink	pink

Verband- und Erste-Hilfe-Scheren  
Bandage- and First Aid Scissors  
145 mm, 5 3/4"



1/1

2/3

BC881R-BC987R

BC881R	schwarz	black	BC983R	gelb	yellow
BC980R	rot	red	BC985R	orange	orange
BC981R	blau	blue	BC987R	pink	pink


Verband- und Erste-Hilfe-Scheren  
Bandage- and First Aid Scissors  
180 mm, 7"



## Gips-Instrumente – Plaster Instruments

Gipsscheren  
Plaster Scissors

X1



MILLIMETERS



**LX548R**

für kleine Gipsverbände  
for small plaster dressings  
160 mm, 6 1/4"



**ESMARCH**

LX550R	200 mm	8"
LX552R	230 mm	9"



**BERGMANN**


**LX553R**  
230 mm, 9"

## Gips-Instrumente – Plaster Instruments

Gipsscheren

Plaster Scissors

LX



BRUNS

LX556R

235 mm, 9 1/4"



LX558R

auch als Kleiderschere verwendbar  
also to be used as bandage scissors  
235 mm, 9 1/4"



SEUTIN

LX565R

auch als Kleiderschere verwendbar  
also to be used as bandage scissors  
230 mm, 9"



## Gips-Instrumente – Plaster Instruments

### Gipsscheren Plaster Shears

#### Besondere Vorteile:

- Die Konstruktion des offenen Kastens verhindert das lästige Verstopfen durch Gipsreste, außerdem ist eine leichte Reinigung der Backen möglich.
- Die tiefe Zahnung an den Backen verhindert das Abgleiten der Gipsschere vom Gips.
- Die Messer, mit Hohlschliff, weisen eine bessere Schnitffähigkeit auf, daher geringerer Kraftaufwand.
- Messerführung durch Kugellager.

#### Special features:

- The construction of the open jaws avoids tiresome clogging by plaster remnants and allows furthermore an easy cleaning.
- The shears do not glide off the plaster owing to the deep teething on the jaws.
- The knives with hollow grinding have a better cutting ability and therefore less force is required.
- The knives are guided by ball bearings.

Kastenteil offen

open jaws



STILLE-AESFULAP

LX574R

260 mm, 10 1/4"




STILLE-AESFULAP

LX575R

360 mm, 14 3/8"

X1



MILLIMETERS

## Gips-Instrumente – Plaster Instruments

Gipsscheren

Plaster Shears

LX



STILLE

LX582R

230 mm, 9"



STILLE

LX587R

Messer mit verstärktem Rücken

Blade with reinforced back


370 mm, 14 2/3"



## Gips-Instrumente – Plaster Instruments

Gipsspreizer, Gipsabreißzangen  
Plaster Spreader, Plaster Breakers

X7



HENNIG

LX600R

270 mm, 10 3/4"



MOD.USA

LX605R

230 mm, 9"



mit abnehmbaren Backen  
with removable jaws

LX610R

290 mm, 11 1/2"



WOLFF

LX618R 180 mm 7"

LX624R 240 mm 9 1/2"

Gipsabreißzangen (Rabenschnabel)  
Plaster breakers

MILLIMETERS